



## Hardness (Total & Calcium) Test Kit, Model HAC-DT

**Product #:** 2063900  
**AUD Price (Incl. GST):** Contact Hach  
**Not in stock - ships within 4-6 weeks**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10 - 4000 mg/L. Contains approximately 100 of each (Total and Calcium) test.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

---

### Specifications

Case Style:	F
Method Name:	Digital Titrator/EDTA
Model:	HAC-DT
Number of tests:	100 Ca +100 T
Parameter:	Hardness, Total & Calcium - as CaCO <sub>3</sub>
Platform :	Digital Titrator
Range:	10 - 4,000 mg/L CaCO <sub>3</sub>
Ship Wt. (lbs):	4
Smallest Increment:	0.1-10

---

### What's in the box?

Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.