



Be Right™



## Sulfite Test Kit, Model SU-DT

**Product #:** 2063300  
AUD Price (Incl. GST): Contact Hach  
Not in stock - ships within 4-6 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Using the Digital Titrator, this kit contains 100 tests with a range of 10-800 mg/L.

**Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Case Style:	F
Method Name:	Digital Titrator/Iodometric
Model:	SU-DT
Number of tests:	100
Parameter:	Sulfite - As $\text{SO}_3^{2-}$
Platform :	Digital Titrator
Range:	4 - 800 mg/L $\text{SO}_3^{2-}$
Ship Wt. (lbs):	4
Smallest Increment Steps:	0.4-4.0**

---

### What's in the box?

Includes Digital Titrator, flask, graduated cylinder, demineralizer bottle, three reagents, instruction manual and carrying case.