



Acidity Test Kit, Model AC-DT

Product #: 2064000
AUD Price (Incl. GST): Contact Hach
Not in stock - ships within 4-6 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Acidity is usually caused by the presence of mineral acids, salts of strong acids, or free carbon dioxide. In surface or ground waters, acidity may form natural substances or industrial pollution. Excessive levels cause corrosion and can be detrimental to fish. Utilizing the Hach Digital Titrator, the Model AC-DT kit is used to determine acidity concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifications

Case Style:	F
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Method Name:	Digital Titrator/Sodium hydroxide
Model:	AC-DT
Number of tests:	100
Parameter:	Acidity - (P or Total) as CaCO_3
Platform :	Digital Titrator
Range:	10 - 4,000 mg/L CaCO_3
Ship Wt. (lbs):	3.5
Smallest Increment Steps:	0.1-1.0

What's in the box?

Includes Digital Titrator, flask, graduated cylinder, four reagents, instruction manual and carrying case.