



Be Right™



Replacement: TUBING, TYGON(R), 3/8" ID X 5/8" OD

Product #: 2943221
AUD Price (Incl. GST): Contact Hach
Not in stock - ships within 4-6 weeks

Replacement Tygon Tubing, 3/8" interior diameter X 5/8" outside diameter for Xenosep SPE Systems Apparatus.

Specifications

Body material:	Tygon
Diameter:	5/8 in
I.D.:	3/8 in
Length:	3 ft
Wall Thickness:	1/4 in