



Be Right™



Beaker, Griffin, Low Form, Glass, 2000 mL

Product #: 50054

AUD Price (Incl. GST): Contact Hach

Not in stock - ships within 4-6 weeks

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifications

Body material:	Glass
Capacity:	2000 mL
Material:	Glass, Griffin, Low Form
Units:	each