



**Be Right™**



## **Beaker, Griffin, Low Form, Glass, 250 mL, 12/pk**

**Product #:** 50076H  
**AUD Price (Incl. GST):** Contact Hach  
**Not in stock - ships within 4-6 weeks**

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### **Specifications**

Body material:	Glass
Capacity:	250 mL
Material:	Glass, Low Form
Units:	12 /pk