



Be Right™



Bottle, Sampling, Polycarbonate, 1000 mL, 4/pk

Product #: 2324663
AUD Price (Incl. GST): Contact Hach
Not in stock - ships within 4-6 weeks

Transparent polycarbonate bottle with cap features wide mouth (3.8cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Specifications

Capacity: 1000 mL
Quantity: 4 /pk
Volume / package size e Post Text: /pk