



Be Right™



ToxTrak™ Toxicity Reagent Set, 25-49 tests

Product #: 2597200
AUD Price (Incl. GST): Contact Hach
Not in stock - ships within 4-6 weeks

ToxTrak™ Toxicity Test Kit

Use the ToxTrak™ Toxicity Reagent Set to analyze your water or wastewater samples for toxicity. Follows Hach Toxicity Method 10017.

The ToxTrak™ reagents perform toxicity screens on drinking water, wastewater and natural water.

Toxicity is compared to a baseline using indigenous biomass or stock bacteria cultures with a range of 0 -100% inhibition.

A colorimetric method is used to determine the toxicity of samples, with analysis done using one of the Hach's trusted spectrophotometers.

Specifications

Description :	Toxicity Reagent set 25-49 tests
Instrument:	All except DR/3000, DR/820, PC II
Method:	10017
Number of tests:	25 - 49
Parameter:	Toxicity
Range:	0 - 100 % Inhibition

What's in the box?

The ToxTrak™ Toxicity Reagent Set includes: Total Bacteria Count media set (15 TBC tubes, 15 individually wrapped sterile pipets), ToxTrak™ powder pillows (pk/50), ToxTrak™ accelerator solution (15mL), disposable culture tubes (5 packs; each pk/6), sterile transfer pipets (pk/15) and sodium thiosulfate standard (100 mL, 0.0246N).

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)